



## Specification Sheet

Rabbit IgG		
Product Name	RB IgG	
Catalog Number	A101-Ab1	
Lot Number		
Description	Rabbit IgG	
Product Profile	Immunogen	N/A
	Specificity	N/A
	Format	Liquid
	Buffer	50 mM NaCl, 20 mM PBS pH 7.4, with 0.1% NaN <sub>3</sub>
	Source/Clone	Normal rabbit sera /polyclonal
	Subclass	
	Purification	Protein A affinity purified from normal rabbit serum
	Concentration	
	Purity	>95%
	Other	
Application	Suitable for use in ELISA or Lateral flow immunoassay. Other applications have not been tested but use in such assays should not be excluded.	
Performance review		
Store at	2-8°C. For long-term storage, freeze at -20°C.	
Precautions and Disclaimer	This product contains NaN <sub>3</sub> , which has been classified as Harmful, in European Directive 67/548/EEC in the concentration range of 0.1–1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent NaN <sub>3</sub> build-up in drains.	
Note	Centrifuge before opening to ensure complete recovery of vial contents	
<b>FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY, NOT FOR USE IN DIAGNOSTIC PROCEDURES</b>		