



Specification Sheet

Goat anti Alpha Human Chorionic Gonadotropin (hCG)		
Product Name		Goat anti-HCG alpha
Catalog Number		A103-Ab1
Lot Number		
Description		polyclonal Goat anti-hCG alpha unit
Product Profile	Immunogen	Human hCG alpha subunit
	Specificity	Specific for human HCG alpha submit
	Format	Liquid
	Buffer	50 mM NaCl, 20 mM PBS pH 7.4, with 0.1% NaN ₃
	Source/Clone	Purified from whole goat anti-sera /polyclonal
	Subclass	N/A
	Purification	Ion-exchange purified
	Concentration	
	Purity	>95%
	Other	
Application		Suitable for use in ELISA or Lateral flow immunoassay as a capture antibody. Other applications have not been tested but use in such assays should not be excluded.
Performance review		<ul style="list-style-type: none"> Great activity
Store at		2-8°C. For long-term storage, freeze at -20°C.
Precautions and Disclaimer		This product contains NaN ₃ , which has been classified as Harmful, in European Directive 67/548/EEC in the concentration range of 0.1–1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent NaN ₃ build-up in drains.
Note		Centrifuge before opening to ensure complete recovery of vial contents
FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY, NOT FOR USE IN DIAGNOSTIC PRODEDURES		