



## Specification Sheet

Antibody against "e" antigen of Hepatitis B Virus		
Product Name	HBeAg mAb 2	
Catalog Number	A28-Ab2	
Lot Number		
Description	Monoclonal antibody to Hepatitis B "e" Antigen (HBeAg)	
Product Profile	Immunogen	Native HBeAg
	Specificity	Specific for HBeAg
	Format	Liquid
	Buffer	50 mM NaCl, 20 mM PBS pH 7.4, with 0.1% NaN <sub>3</sub>
	Source/Clone	Mouse ascites / Monoclonal
	Subclass	IgG1
	Purification	Ion-exchange chromatography purified
	Concentration	
	Purity	>95%
Application	Suitable for use in ELISA or Lateral flow immunoassay for HBeAg and HBeAb detection. Other applications have not been tested but use in such assays should not be excluded.	
Performance review		
Store at	2-8°C. For long-term storage, freeze at -20°C.	
Precautions and Disclaimer	This product contains NaN <sub>3</sub> , which has been classified as Harmful, in European Directive 67/548/EEC in the concentration range of 0.1–1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent NaN <sub>3</sub> build-up in drains.	
Note	Centrifuge before opening to ensure complete recovery of vial contents	
<b>FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY, NOT FOR USE IN DIAGNOSTIC PROCEDURES</b>		