



Specification Sheet

Antibody against <i>Vibrio Cholerae O1-Ogawa</i>		
Product Name	VCO1 Ogawa mAb1	
Catalog Number	A41-Ab1	
Lot Number		
Description	Monoclonal antibody to <i>Vibrio Cholerae O1-Ogawa</i>	
Product Profile	Immunogen	inactive <i>whole Vibrio Cholerae O1-Ogawa</i>
	Specificity	Specific for <i>Vibrio Cholerae O1-Ogawa</i> Negative for related bacteria (<i>Vibrio Cholerae O139,, O1-Inaba, etc</i>)
	Format	Liquid
	Buffer	50 mM NaCl, 20 mM PBS pH 7.4, with 0.1% NaN ₃
	Source/Clone	Mouse ascites / Monoclonal
	Subclass	IgG2a
	Purification	Protein A affinity purified
	Concentration	
Purity	>95%	
Application	Suitable for use in ELISA or Lateral flow immunoassay (double antibody test for <i>Vibrio Cholerae O1-Ogawa</i> detection). Other applications have not been tested but use in such assays should not be excluded.	
Performance review	<ul style="list-style-type: none"> Capable of creating tests with sensitivity at 10⁴cfu/ml or lower Recommended pair for double antibody immunoassay: Coating: A41-Ab2; labeling: A41-Ab1 	
Store at	2-8°C. For long-term storage, freeze at -20°C.	
Precautions and Disclaimer	This product contains NaN ₃ , which has been classified as Harmful, in European Directive 67/548/EEC in the concentration range of 0.1–1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent NaN ₃ build-up in drains.	
Note	Centrifuge before opening to ensure complete recovery of vial contents	
FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY, NOT FOR USE IN DIAGNOSTIC PRODEDURES		