



Specification Sheet

Antibody against <i>Vibrio Cholerae O1 Inaba</i>		
Product Name	VCO1 Inaba mAb2	
Catalog Number	A42-Ab2	
Quantity		
Lot Number		
Description	Monoclonal antibody to <i>Vibrio Cholerae O1 Inaba</i>	
Product Profile	Immunogen	inactive whole <i>Vibrio Cholerae O1-Inaba</i>
	Specificity	Specific for <i>Vibrio Cholerae O1- Inaba</i> Cross with <i>Vibrio Cholerae O1- Ogawa</i> Negative for related bacteria (<i>Vibrio Cholerae O139, etc</i>)
	Format	Liquid
	Buffer	50 mM NaCl, 20 mM PBS pH 7.4, with 0.1% NaN ₃
	Source/Clone	Mouse ascites / Monoclonal
	Subclass	IgG2b
	Purification	Protein A affinity purified
	Concentration	
Purity	>95%	
Application	Suitable for use in ELISA or Lateral flow immunoassay (double antibody test for <i>Vibrio Cholerae O1-inaba</i> detection). Other applications have not been tested but use in such assays should not be excluded.	
Performance review	<ul style="list-style-type: none"> Capable of creating tests with sensitivity at 10³cfu/ml or lower Recommended pair for double antibody immunoassay: Coating: A42-Ab2; labeling: A42-Ab1 	
Store at	2-8°C. For long-term storage, freeze at -20°C.	
Precautions and Disclaimer	This product contains NaN ₃ , which has been classified as Harmful, in European Directive 67/548/EEC in the concentration range of 0.1–1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent NaN ₃ build-up in drains.	
Note	Centrifuge before opening to ensure complete recovery of vial contents	
FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY, NOT FOR USE IN DIAGNOSTIC PRODEDURES		