



## Specification Sheet

Antibody against <i>Vibrio Cholerae O1</i>		
Product Name	VCO1 mAb3	
Catalog Number	A43-Ab3	
Lot Number		
Description	Monoclonal antibody to <i>Vibrio Cholerae O1</i>	
Product Profile	Immunogen	inactive whole <i>Vibrio Cholerae O1-Inaba</i>
	Specificity	Specific for <i>Vibrio Cholerae O1-Ogawa and Inaba</i> Negative for related bacteria ( <i>Vibrio Cholerae O139, etc</i> )
	Format	Liquid
	Buffer	50 mM NaCl, 20 mM PBS pH 7.4, with 0.1% NaN <sub>3</sub>
	Source/Clone	Mouse ascites / Monoclonal
	Subclass	IgG1
	Purification	Protein A affinity purified
	Concentration	
Purity	>95%	
Application	Suitable for use in ELISA or Lateral flow immunoassay (double antibody test for <i>Vibrio Cholerae O1-Ogawa</i> detection). Other applications have not been tested but use in such assays should not be excluded.	
Performance review		
Store at	2-8°C. For long-term storage, freeze at -20°C.	
Precautions and Disclaimer	This product contains NaN <sub>3</sub> , which has been classified as Harmful, in European Directive 67/548/EEC in the concentration range of 0.1–1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent NaN <sub>3</sub> build-up in drains.	
Note	Centrifuge before opening to ensure complete recovery of vial contents	
<b>FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY, NOT FOR USE IN DIAGNOSTIC PROCEDURES</b>		