



Specification Sheet

Antibody against Dengue NS1 (Type I-IV)		
Product Name		NS1 mAb 6
Catalog Number		A47-Ab6
Quantity		
Lot Number		
Description		Monoclonal antibody to Dengue NS1 Protein (Type I-IV)
Product Profile	Specificity	Specific for Dengue NS1, Type I-IV
	Immunogen	Native Dengue NS1 Antigen
	Format	Liquid
	Buffer	50 mM NaCl, 20 mM PBS pH 7.4, with 0.1% NaN ₃
	Source/Clone	Mouse ascites / Monoclonal
	Subclass	IgG1
	Purification	Affinity Purification
	Concentration	
	Purity	>90%
	Affinity constant	
Application		Suitable for labeling in ELISA or Lateral flow immunoassay (double antibody test for Dengue NS1 antigen detection). Other applications have not been tested but use in such assays should not be excluded.
Performance review		<ul style="list-style-type: none">• Suitable also for western Blot• Recommended pair for double antibody immunoassay: Coating: A51-Ab6; labeling: A51-Ab5
Store at		Store at: 2-8°C for less than six months, -20°C for over six months. Avoid multiple freeze/thaw cycles.
Precautions and Disclaimer		This product contains NaN ₃ , which has been classified as Harmful, in European Directive 67/548/EEC in the concentration range of 0.1–1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent NaN ₃ build-up in drains.
Note		Centrifuge before opening to ensure complete recovery of vial contents
FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY, NOT FOR USE IN DIAGNOSTIC PRODEDURES		

Artron BioResearch Inc
3938 North Fraser Way
Burnaby, BC Canada V5J 5H6

Tel: +1 (604) 415-9757
Fax: +1 (604) 415-9795

info@artronbio.com
www.artronbio.com