



Specification Sheet

Rabbit anti Streptococcus A		
Product Name	STA mAb 1	
Catalog Number	A50-Ab1	
Quantity		
Lot Number		
Description	Rabbit polyclonal antibody to Streptococcus A	
Product Profile	Specificity	Specific for Streptococcus A
	Immunogen	Intact molecule Streptococcus A
	Format	Liquid
	Buffer	50 mM NaCl, 20 mM PBS pH 7.4, with 0.1% NaN ₃
	Source/Clone	Rabbit / Polyclonal
	Subclass	N/A
	Purification	Ion-exchange purified
	Concentration	
	Purity	>95%
Application	Suitable for use in ELISA or Lateral flow immunoassay (double antibody test for AFP detection). Other applications have not been tested but use in such assays should not be excluded.	
Performance review	<ul style="list-style-type: none"> Specificity tested by using a complete panel of organisms obtained from the ATCC. Recommended pair for double antibody immunoassay: Coating: A51-Ab2; labeling: A51-Ab1 	
Store at	2-8°C. For long-term storage, freeze at -20°C.	
Precautions and Disclaimer	This product contains NaN ₃ , which has been classified as Harmful, in European Directive 67/548/EEC in the concentration range of 0.1–1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent NaN ₃ build-up in drains.	
Note	Centrifuge before opening to ensure complete recovery of vial contents	
FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY, NOT FOR USE IN DIAGNOSTIC PROCEDURES		