



## Specification Sheet

Antibody against alpha -Fetoprotein		
Product Name	AFP mAb2- coating	
Catalog Number	A51-Ab2	
Lot Number		
Description	Monoclonal antibody to Human alpha -Fetoprotein	
Product Profile	Specificity	Specific for alpha -Fetoprotein
	Immunogen	Purified native Human AFP protein
	Format	Liquid
	Buffer	50 mM NaCl, 20 mM PBS pH 7.4, with 0.1% NaN <sub>3</sub>
	Source/Clone	Mouse ascites / Monoclonal
	Subclass	IgG1
	Purification	Ion-exchange purified
	Concentration	
	Purity	>95%
Application	Suitable for use in ELISA or Lateral flow immunoassay (double antibody test for AFP detection). Other applications have not been tested but use in such assays should not be excluded.	
Performance review	<ul style="list-style-type: none"> <li>Capable of creating tests with sensitivity of 25 ng/ml or lower.</li> <li>Recommended pair for double antibody immunoassay: Coating: A51-Ab2; labeling: A51-Ab1</li> </ul>	
Store at	2-8°C. For long-term storage, freeze at -20°C.	
Precautions and Disclaimer	This product contains NaN <sub>3</sub> , which has been classified as Harmful, in European Directive 67/548/EEC in the concentration range of 0.1–1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent NaN <sub>3</sub> build-up in drains.	
Note	Centrifuge before opening to ensure complete recovery of vial contents	
<b>FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY, NOT FOR USE IN DIAGNOSTIC PRODEDURES</b>		