



Specification Sheet

Antibody against Human Carcinoembryonic Antigen		
Product Name	CEA mAb1-Labeling	
Catalog Number	A52-Ab1	
Lot Number		
Description	Monoclonal antibody to Human Carcinoembryonic Antigen	
Product Profile	Specificity	Specific for Human Carcinoembryonic Antigen Negative for human albumin, NCA, NCA-2, and PSA, etc
	Format	Liquid
	Buffer	50 mM NaCl, 20 mM PBS pH 7.4, with 0.1% NaN ₃
	Source/Clone	Mouse ascites / Monoclonal
	Subclass	IgG1
	Purification	Ion-exchange purified
	Concentration	
	Purity	>95%
	Immunogen	CEA purified from CEA culture supernatant
Application	Suitable for use in ELISA or Lateral flow immunoassay (double antibody test for CEA detection). Other applications have not been tested but use in such assays should not be excluded.	
Performance review	<ul style="list-style-type: none"> Capable of creating tests with sensitivity of 5 ng/ml or lower. Recommended pair for double antibody immunoassay: Coating: A52-Ab2; labeling: A52-Ab1 	
Store at	2-8°C. For long-term storage, freeze at -20°C.	
Precautions and Disclaimer	This product contains NaN ₃ , which has been classified as Harmful, in European Directive 67/548/EEC in the concentration range of 0.1–1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent NaN ₃ build-up in drains.	
Note	Centrifuge before opening to ensure complete recovery of vial contents	
FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY, NOT FOR USE IN DIAGNOSTIC PROCEDURES		