



Specification Sheet

Antibody against Human hemoglobin		
Product Name	FOB mAb3-labeling	
Catalog Number	A54-Ab3	
Lot Number		
Description	Monoclonal antibody to human hemoglobin	
Product Profile	Specificity	Specific for Human hemoglobin
	Format	Liquid
	Buffer	50 mM NaCl, 20 mM PBS pH 7.4, with 0.1% NaN ₃
	Source/Clone	Mouse ascites / Monoclonal
	Subclass	IgG1
	Purification	Ion-exchange purified
	Concentration	
	Purity	>95%
Affinity constant		
Application	Suitable for use in ELISA or Lateral flow immunoassay (double antibody test for FOB detection). Other applications have not been tested but use in such assays should not be excluded.	
Performance review	<ul style="list-style-type: none"> Capable of creating tests with sensitivity of 30 ng/ml or lower. No cross reaction with hemoglobin from other species. Recommended pair for double antibody immunoassay: Coating: A54-Ab4; labeling: A54-Ab3 	
Store at	2-8°C. For long-term storage, freeze at -20°C.	
Precautions and Disclaimer	This product contains NaN ₃ , which has been classified as Harmful, in European Directive 67/548/EEC in the concentration range of 0.1–1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent NaN ₃ build-up in drains.	
Note	Centrifuge before opening to ensure complete recovery of vial contents	
FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY, NOT FOR USE IN DIAGNOSTIC PRODEDURES		