



Specification Sheet

Antibody against Human Myoglobin	
Product Name	Myo mAb2-coating
Catalog Number	A61-Ab2
Lot Number	
Description	Monoclonal antibody to human myoglobin
Specificity	Specific for Human myoglobin
Format	Liquid
Buffer	50 mM NaCl, 20 mM PBS pH 7.4, with 0.1% NaN ₃
Source/Clone	Mouse ascites / Monoclonal
Subclass	IgG2b
Purification	Protein A affinity purified
Concentration	
Purity	>95%
other	
Application	Suitable for use in ELISA or Lateral flow immunoassay (double antibody test for myoglobin detection). Other applications have not been tested but use in such assays should not be excluded.
Performance review	<ul style="list-style-type: none"> Recommended pair for double antibody immunoassay: Coating: A61-Ab2; labeling: A61-Ab1
Store at	2-8°C. For long-term storage, freeze at -20°C.
Precautions and Disclaimer	This product contains NaN ₃ , which has been classified as Harmful, in European Directive 67/548/EEC in the concentration range of 0.1–1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent NaN ₃ build-up in drains.
Note	Centrifuge before opening to ensure complete recovery of vial contents
FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY, NOT FOR USE IN DIAGNOSTIC PROCEDURES	