



## Specification Sheet

Antibody against Human Fatty Acid Binding Protein		
Product Name	FABP mAb1-Labeling	
Catalog Number	A63-Ab1	
Lot Number		
Description	Monoclonal antibody to human fatty acid binding protein	
Product Profile	Immunogen	Purified human FABP3
	Specificity	Specific for human fatty acid binding protein
	Format	Liquid
	Buffer	50 mM NaCl, 20 mM PBS pH 7.4, with 0.1% NaN <sub>3</sub>
	Source/Clone	Mouse ascites / Monoclonal
	Subclass	IgG1
	Purification	Protein A affinity purified
	Concentration	
	Purity	>95%
Affinity constant		
Application	Suitable for use in ELISA or Lateral flow immunoassay (double antibody test for FABP detection). Other applications have not been tested but use in such assays should not be excluded.	
Performance review	<ul style="list-style-type: none"> <li>Recommended pair for double antibody immunoassay: Coating: A63-Ab2; labeling: A63-Ab1</li> </ul>	
Store at	2-8°C. For long-term storage, freeze at -20°C.	
Precautions and Disclaimer	This product contains NaN <sub>3</sub> , which has been classified as Harmful, in European Directive 67/548/EEC in the concentration range of 0.1–1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent NaN <sub>3</sub> build-up in drains.	
Note	Centrifuge before opening to ensure complete recovery of vial contents	
<b>FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY, NOT FOR USE IN DIAGNOSTIC PROCEDURES</b>		