



Specification Sheet

Antibody against Human IgA		
Product Name	Anti-Human IgA mAb2	
Catalog Number	A92-Ab2	
Lot Number		
Description	Monoclonal antibody to Human IgA	
Product Profile	Immunogen	Purified Human Serum IgA
	Specificity	Specific for Human IgA No cross reactivity with human IgM, IgD, IgE or IgG.
	Format	Liquid
	Buffer	50 mM NaCl, 20 mM PBS pH 7.4, with 0.1% NaN ₃
	Source/Clone	Mouse ascites / Monoclonal
	Subclass	IgG1
	Purification	Ion-exchange purified
	Concentration	
	Purity	>95%
Application	Suitable for use in Western Blot, ELISA or Lateral flow immunoassay as a capture antibody. Other applications have not been tested but use in such assays should not be excluded.	
Performance review	N/A	
Store at	2-8°C. For long-term storage, freeze at -20°C.	
Precautions and Disclaimer	This product contains NaN ₃ , which has been classified as Harmful, in European Directive 67/548/EEC in the concentration range of 0.1–1.0%. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent NaN ₃ build-up in drains.	
Note	Centrifuge before opening to ensure complete recovery of vial contents	
FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY, NOT FOR USE IN DIAGNOSTIC PROCEDURES		