



Certificate of Analysis

Antibody against Human red blood cell		
Product Name	RBC ABO mAb1	
Catalog Number	A96-Ab1	
Lot Number		
Description	Monoclonal antibody to human red blood cell (ABO)	
Product Profile	Specificity	Specific for Human red blood cell (ABO)
	Format	Liquid
	Buffer	50 mM NaCl, 20 mM PBS pH 7.4, with 0.1% NaN ₃
	Source/Clone	Mouse ascites / Monoclonal
	Subclass	to be determined
	Purification	Ion-exchange chromatography purified
	Concentration	
	Purity	>90%
Application	Suitable for removing blood red cell in whole blood with usage 0.5-1mg/ml (final concentration) in sample pad.	
Store at	2-8°C for less than six months, -20°C for over six months. Avoid multiple freeze/thaw cycles.	
Note	Centrifuge before opening to ensure complete recovery of vial contents	
FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY, NOT FOR USE IN DIAGNOSTIC PROCEDURES		